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**Intended use**

- conservative treatment in motor segment dysfunctions of different aetiology and severity.
- stabilisation of the lumbar region in cases of segmental instability following surgery.
- recovery from neurosurgery and orthopaedic surgery – functional treatment.
- prevention of re-injury, occupational and sport overloads.
- congenital defects of the lumbar spine (with or without pain syndromes).
- chronic pain in the course of degenerative disease.

Indications:

- supporting the treatment and rehabilitation process in the case of injuries, diseases and dysfunctions of the intervertebral joints, intervertebral disc, muscle attachments and ligaments of the thoracolumbar section and lumbosacral transition of the spine.
- controlling range of motion and correcting posture

Contraindications:

- hard-to-heal wounds or skin lesions within the area of the orthosis;
- allergic reaction to the materials used;

Target groups of patients:

The orthosis is intended for patients after injury and surgery. It supports the treatment and rehabilitation process, leading to the alleviation of clinical symptoms associated with degenerative changes and functional disorders of the above-mentioned motor segment.

Disposal of the medical device:

The device must be disposed of following the applicable local environmental regulations.

Choosing your orthosis:

1. The first time the orthosis is fitted, it should be done upon prescription and under the supervision of qualified medical personnel (doctor or physiotherapist).
2. The orthosis should be selected based on the dimensions shown in the size chart.
3. If the device becomes too loose or changes its position, it must be refitted starting from the first step in the section „How to put the orthosis on”.

Size chart:

Size	1	2	3	4	5	6
Waist circumference (cm)	52 - 62	63 - 74	75 - 89	90 - 105	106 - 120	121 - 140

How to put the orthosis on:

Optimal fitting of the orthosis by qualified medical personnel is essential to ensure proper stabilisation. The orthosis should be put on in the supine position to allow for easy and proper adjustment.

1. Put the orthosis on so that its bottom edge is about 2 cm above the end of the tailbone (two fingers' width from the gluteal cleft) and the tag sewn on the inside is on top. The back reinforcements (liners) should cling to the thoracic-lumbosacral region.
2. Fasten the tension strap first and then fasten the tightening straps.
3. Use the tightening strap to adjust the dorsal pressure. The orthosis should be fastened with a force that provides stabilisation (partially restricts spinal movement in the lumbosacral segment, but does not restrict free breathing).
4. The duration of use of the orthosis should be strictly prescribed by a doctor or physiotherapist. A rehabilitation exercise programme should be introduced in parallel.

Instructions for cleaning:

The orthosis must be washed by hand at low temperatures (30° C) with mild detergents, without softeners. Do not use bleach or strong detergents. Do not spin, do not wring, gently squeeze before drying. To dry, spread the parts of the upholstery on a towel and air-dry. Keep away from artificial sources of heat. Check the condition of the orthosis after each wash. Clean the product systematically, at least twice a month, and more often if necessary.

Composition:

PA, PES, PE

Intended user:

The medical device is intended for regular (non-professional user) as well as for the professional user.

**ATTENTION!**

Do not use the device if you have circulatory disorders.
Do not use the device if it has been damaged or deformed.
Do not attempt to repair the device yourself.
Do not use the device of the wrong size or in a way for which it is not intended. If you experience any problems related to using the product, please get in touch with your doctor. Failure to observe proper personal hygiene and inappropriate laundry detergents can cause abrasions or other ailments associated with increased sweating and bacterial growth. Do not use the orthosis with warming gels or ointments. Any serious incident in connection with a device must be reported to the manufacturer and the competent authority of the Member State where the user and/or patient is based.

IT IS RECOMMENDED TO CONSULT A SPECIALIST BEFORE USING THE DEVICE

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